TaKaRa EpiTaq™ HS (for bisulfite-treated DNA)

Code No. R110A
Size : 250 units

Shipping at — 20°C
Store at — 20°C

Supplied Reagents :

10X EpiTaq PCR Buffer (Mg²⁺ free) 1 ml
dNTP Mixture (2.5 mM each) 1.2 ml
25 mM MgCl₂ 1.2 ml

Lot No.

Concentration : 5 units/μl

Volume : 50 μl

Expiration Date:

For a detailed protocol, refer to the product user manual.

Storage Buffer : 20 mM Tris-HCl (pH8.0)
100 mM KCl
0.1 mM EDTA
1 mM DTT
0.5% Tween 20
0.5% Nonidet P-40
50% Glycerol

Supplied dNTP Mixture (2.5 mM each)

dNTP Mixture is ready for use in PCR without dilution.

Form : Dissolved in water (sodium salts), pH7 - 9

Purity : ≥ 98% for each dNTP

Unit definition : One unit is the amount of enzyme that will incorporate 10 nmol of dNTPs into acid-insoluble products in 30 min. at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition :

25 mM TAPS (pH9.3 at 25°C)
50 mM KCl
2 mM MgCl₂
0.1 mM DTT
200 μM each dATP·dGTP·dCTP
100 μM [³H]-dTTP
0.25 mg/ml activated salmon sperm DNA

Purity : Nicking, endonuclease, and exonuclease activity were not detected after incubation of 0.6 μg of supercoiled pBR322 DNA, 0.6 μg of λ DNA, or 0.6 μg of λ-Hind III digest with 10 units of this enzyme for 1 hour at 74°C.

Application : PCR amplification using bisulfite-treated DNA as a template.

PCR products : As most PCR products amplified with TaKaRa EpiTaq HS have one A at the 3'-termini, the obtained PCR products can be directly cloned into T-vectors. Also it is possible to clone the products in blunt end vectors after blunting and phosphorylation of the ends.

PCR test :

Good performance by PCR was confirmed by amplification of a 148-bp fragment from a DNA template containing dUTP. Antibody-mediated inhibition of EpiTaq polymerase activity was confirmed to be more than 90% after the reaction was incubated at 55°C for 10 min.

General reaction mixture for PCR (total 50 μl) :

	Volume	Final conc.
TaKaRa EpiTaq HS (5 U/μl)	0.25 μl	1.25 U / 50 μl
10X EpiTaq PCR Buffer (Mg ²⁺ free)	5 μl	1X
25 mM MgCl ₂	5 μl	2.5 mM
dNTP Mixture (2.5 mM each)	6 μl	0.3 mM
Template	< 100 ng	
Primer 1	20 pmol	0.4 μM
Primer 2	20 pmol	0.4 μM
Sterile distilled water	up to 50 μl	

Try the above reaction mixture first. If necessary, change the concentration of MgCl₂, dNTP Mixture, or primers to improve amplification efficiency, extension, or specificity. Please refer to the "Troubleshooting" section in the product user manual for additional details.

PCR conditions : This enzyme can be used with standard PCR conditions, since the monoclonal antibody is denatured in the initial DNA denaturation step. There is no need for an additional step to denature the antibody to Taq polymerase.

Example : Amplification of a 0.5 kb DNA fragment

98°C 10 sec. *1
55°C 30 sec.
72°C 30 sec. for <500 bp or 1 min. for 500 bp - 1 kb *2

30 - 40 cycles

*1 : Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. Denaturation for 5 - 10 sec. at 98°C, or 20 - 30 sec. at 94°C is recommended.

*2 : Please set it according to the amplicon size. The standard setting is 30 sec. for products < 500 bp and 1 min. for products 500 bp - 1 kb in size. For amplification of targets longer than 1 kb, the standard setting is 1 min./kb.

NOTICE TO PURCHASER: LIMITED LICENSE

[P1] PCR Notice

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,789,224, 5,618,711, 6,127,155 and claims outside the US corresponding to expired US Patent No. 5,079,352. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

EpiTaq is a trademark of TAKARA BIO INC.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

TaKaRa EpiTaq™ HS (for bisulfite-treated DNA)

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